



CS Ceramic CO.,Ltd

Tel.: 0086-731-23455639 Web: www.csceramic.com

Contact: Ms.Eve Zhang (eve@csceramic.com)

Name: carbon sulfur crucible / Ceramic crucibles for LECO analysis instrument

Size: 25 X 25mm

Blank value of Carbon and sulfur element : $C \leq 3\text{ppm}$, $S \leq 2\text{ppm}$

Application: test carbon and sulphur elements of iron and steel material

MOQ: 10 cartons, a carton has 1,000 pieces crucibles

Weight : 20KG /carton;1000pcs/carton ;L 53.5 X W26.5 X H16.5 cm

Payment Terms: TT 100% before loading

Delivery time: stock

Packing :normal packing+plywood packing



Packing photos:



Compatible Instruments:

Eltra CS2000

Eltra CS800

Eltra Helios

Horiba EMIA-320V2/920V2

Horiba EMIA820V/320V/220V

Leco CS125/225

Leco CS200/C200/S200/WC200

Leco CS230

Leco CS244/344/IR212/232

Leco CS300

Leco CS400/C400/S400

Leco CS444/IR412/432

Leco CS600

Leco WR112

OEM Equivalent Parts:

Alpha Resources: AR 3818

Eltra: 90150

Horiba: 905.200.380.001AN

Horiba: 905.200.380.001

Leco: 528-018

SerCon: SC0893

Direction for use:

To get a more exact test result, before use it, suggest baking the crucibles under 1100 degrees centigrade for 2 hours, then the blank of C&S will be very stable and keep low PPM.

Company introduction:

CS Ceramic CO.,Ltd is a manufacturer working in the field since 1974, it is the standard maker for carbon-sulfur crucibles in China. Other consumables for LECO are also supplied. OEM is welcome as we have R&D department.

Any questions or ideas, please feel free to contact me.

Ms.Eve Zhang (eve@csceramic.com)

Tel: 0086-731-23455639

Website: www.csceramic.com



Linkdin: <http://www.linkedin.com/company/cs-ceramic-co-ltd?trk=biz-companies-cym>